


**Sample Name: Hydroxy Terminated  
Polyisoprene, 1,4-microstructure**

**Sample #: P6202-IPOH**

**Structure:**




**Composition:**

Mn x 10 <sup>3</sup>	PDI
26.0	1.05

**Synthesis Procedure:**

1,4-addition hydroxy terminated polyisoprene was prepared by anionic living polymerization in a non-polar solvent followed by termination with ethylene oxide. The scheme of the reaction is illustrated below:



**Characterization:**

The molecular weight and polydispersity index (PDI) are obtained by size exclusion chromatography (SEC) in THF. SEC analysis was performed on a Varian liquid chromatograph equipped with refractive and UV light scattering detectors. Three SEC columns from Supelco (G6000-4000-2000 HXL) were used with triple detectors from Viscotek Co.

**Functionality:** functionality of the obtained polymer was determined by reacting polymer in dried non quantity of acetic anhydride in the presence of pyridine as a catalyst and the liberated COOH was titrated by acid-base titration.

**Solubility:**

Hydroxy terminated polyisoprene is soluble in DMF, THF, toluene, hexane, cyclohexane and CHCl<sub>3</sub>. It precipitates from methanol, ethanol, water.

**SEC of Sample:**

