


**Sample Name:** Poly(2-chloro styrene)

**SEC of Homopolymer:**

**Sample #:** P6661B-2CIS

**Structure:**




**Composition:**

$M_n \times 10^3$	PDI
30.4	1.65

**Synthesis Procedure:**

Poly(2-chloro styrene) is synthesized by controlled radical polymerization of 2-chloro styrene and the reaction scheme is shown below.



**Characterization:**

The molecular weight and polydispersity index (PDI) are obtained by size exclusion chromatography (SEC) in THF. SEC analysis was performed on a Varian liquid chromatograph equipped with refractive and UV light scattering detectors. Three SEC columns from Supelco (G6000-4000-2000 HXL) were used with triple detectors from Viscotek Co.

**Solubility:**

Poly(2-chloro styrene) is soluble in DMF, THF, toluene and  $\text{CHCl}_3$ . It precipitates from methanol, ethanol, water and hexanes.

